**Add this addendum to the front of your HAMILTON-C2 Operator's Manual.**

The HAMILTON-C2 software version 2.2.x introduces some important enhancements and updates to the device software to further improve the performance of the device and changes to the Operator's Manual. This addendum summarizes the changes.

**Changes to Operator's Manual**

This addendum provides updates for the *HAMILTON-C2 Operator's Manual*, PN 624131/07, for software version 2.2.x. It reflects the new functionality and options, and includes updated safety information. Review the information before operating the device with the updated software.

The following sections are organized by chapters in the *HAMILTON-C2 Operator's Manual* with related information for ease of reference.

**What's new in software version 2.2.x and the operator's manual**

The following features and options have been added or updated in version 2.2.x.

**Table 1-1. Changes in software version 2.2.x and the operator's manual**

<b>Feature/Option</b>	<b>Additional information</b>
New display languages added: Croatian, Finnish, Romanian	--
Once activated, CO2 monitoring remains activated until manually deactivated	--
Changed the method for calibrating the infant flow sensor; a calibration adapter is now used	See "Chapter 3. Tests, calibrations, and utilities" on page 3 of this addendum
Alarm loudness setting moved to Settings window and Alarm limit settings and ranges table changed	See "Chapter 4. Ventilator settings" on page 4 of this addendum
Blower (turbine) lifetime updated	Chapter 10, Maintenance The service interval/lifetime for the turbine has been extended as described here. Interval: Dynamic lifetime surveillance. Typically 8 years.
Default values changed for the following parameters: ETS, I:E, Rate, TI, Vt/kg	See "Appendix A. Specifications" on page 8 of this addendum

**Table 1-1. Changes in software version 2.2.x and the operator's manual**

<b>Feature/Option</b>	<b>Additional information</b>
P-ramp limitation by one-third TI added	See "Appendix A. Specifications" on page 8 of this addendum
HPO/LPO selection moved to Utilities window, available in Standby	See "Appendix J. Configuration" on page 13 of this addendum
<p>The following alarms have been added or updated:</p> <ul style="list-style-type: none"> <li>• Battery low</li> <li>• Loss of external power</li> <li>• Loss of PEEP</li> <li>• Performance limited by altitude</li> <li>• Battery replacement required</li> <li>• High PEEP</li> <li>• High oxygen</li> <li>• Inspiratory volume limitation (now disabled in noninvasive modes)</li> <li>• CO2 sensor warmup</li> <li>• Function key not operational</li> <li>• Invalid option board</li> <li>• Replace O2 cell</li> <li>• Suctioning maneuver</li> <li>• Technical state failed</li> <li>• Touch not functional</li> <li>• Unknown part number</li> <li>• Release valve defective</li> </ul>	See "Chapter 8. Responding to alarms" on page 4 of this addendum
<p>Replaced the numerical code with text and added troubleshooting information to the Operator's Manual for the following technical alarms:</p> <ul style="list-style-type: none"> <li>• Battery communication error (244011)</li> <li>• Battery totally discharged (444001)</li> <li>• Blower fault (431001)</li> <li>• Buzzer defective (243004)</li> <li>• Loudspeaker defective (243003)</li> <li>• Options not found (249001)</li> <li>• Self test failed (485002)</li> </ul>	See "Chapter 8. Responding to alarms" on page 4 of this addendum
Ability to configure a minimum alarm volume (loudness) for the device	See "Setting the minimum alarm loudness" on page 13 of this addendum

## Chapter 3. Tests, calibrations, and utilities

The calibration method for the infant flow sensor has changed. Note that this applies only to neonatal ventilation.


### Calibrating the infant flow sensor

#### NOTE


An infant flow sensor is required with breathing circuits used for all ventilation modes.

#### To calibrate the infant flow sensor


- 1 Set the ventilator up as for normal ventilation, complete with breathing circuit and expiratory membrane and cover.
- 2 Make sure that the Neonatal patient group is selected, an infant flow sensor is installed, and the calibration adapter is available.
- 3 In the System -> Tests & calib window, select **Flow Sensor**.  
If you have not already disconnected the patient, the text **Disconnect patient** is displayed.
- 4 Disconnect the patient now.


- 5 Follow the instructions displayed in the message line:
  - a. Attach the calibration adapter to the patient end of the flow sensor.

- b. When prompted, turn the flow sensor around as indicated and attach the calibration end to the Y-piece.

- c. When prompted to turn the flow sensor again, turn the flow sensor back to its starting position, and remove the calibration adapter.

- 6 When calibration is complete, verify that there is a green check mark in the **Flow Sensor** checkbox.
- 7 If the calibration is successful, connect the patient, and touch the **Start ventilation** button in the Standby window to start ventilation.

## Chapter 4. Ventilator settings

### Adjusting the alarm loudness

---

**WARNING**

Be sure to set the auditory alarm volume above the ambient sound level. Failure to do so can prevent you from hearing and recognizing alarm conditions.

---

**NOTE**

- The alarm volume cannot be set lower than the minimum alarm loudness specified for the device in Configuration.
  - If the alarm loudness was set to a value below the default setting (5 for adult/pediatric, 3 for neonates) before the ventilator was turned off, it will be reset to the default value when the ventilator is turned back on.
  - If you decrease the alarm loudness during the night shift, do not forget to return it to its daytime setting!
  - The alarm volume control is on the **Settings** tab.
- 

**To adjust the alarm volume**

- 1 Open the **System** -> **Settings** -> **Loudness** window.
- 2 Adjust the **Loudness** value as desired.
- 3 **Test** the loudness as desired.

## Chapter 8. Responding to alarms

Table 8-2 has been updated. The following table lists only those alarms that have been updated or added.

Table 8-2. Alarms and other messages

Alarm	Definition	Action needed
Battery communication error	<i>High priority.</i> Battery data is not available. Ventilation continues.	Make sure the battery connectors are intact and the battery is installed correctly.  If the problem persists, replace the battery.  If the problem still persists, have the ventilator serviced.
Battery low	The low battery alarm has different levels of priority, depending on how much charge is left, and which power supply is in use.  Note that at 20% battery charge, the ventilator can generally continue operation for up to approx 10 min, depending on battery and operating conditions.  <i>High priority.</i> The ventilator is running on battery power, and the total battery charge is below 15%.  <i>Medium priority.</i> The ventilator is running on battery power, and the total battery charge is below 20%.  <i>Low priority.</i> The ventilator is running on AC/DC power, and the total battery charge is below 20%.	Connect the ventilator to its primary power source.  Install charged batteries.
Battery 1+ 2: replacement required	<i>Low priority.</i> Battery capacity is insufficient for reliable operation and must be replaced immediately.	Replace the battery.
<b>NOTE</b> Battery life indications are approximate. The actual battery life depends on ventilator settings, battery age, and level of battery charge.		
Battery totally discharged	<i>High priority.</i> The battery charge level is below 5%.  The ventilator switches to the ambient state.	Connect the device to the primary power supply and recharge the battery.  Provide alternative ventilation.  Have the ventilator serviced.

**Table 8-2. Alarms and other messages**

<b>Alarm</b>	<b>Definition</b>	<b>Action needed</b>
Blower fault	<i>High priority.</i> A blower malfunction was detected. A technical alarm cannot typically be corrected by the operator. The ventilator switches to the ambient state.	Provide alternative ventilation. Have the ventilator serviced.
Buzzer defective	<i>High priority.</i> A buzzer malfunction was detected. A technical alarm cannot typically be corrected by the operator.	Restart device. If the problem persists, have the ventilator serviced.
CO2 sensor warmup	<i>Low priority.</i> CO2 operating temperature not yet reached or unstable.	Wait for sensor to warm up
Function key not operational	<i>Medium priority.</i> Function key defective.	Have the ventilator serviced.
High oxygen	<i>High priority.</i> With LPO selected: The measured oxygen is > the set high Oxygen alarm limit. With HPO selected: The measured oxygen is > 5% over the Oxygen control setting.	Calibrate the oxygen cell. Install a new oxygen cell. Check alarm limits.
High PEEP	<i>Medium priority.</i> Monitored PEEP > (set PEEP + 5 cmH2O) for two consecutive breaths. For DuoPAP and APRV only: <ul style="list-style-type: none"> <li>Alarm applies to both Phigh and Plow settings. The alarm sounds when the monitored Phigh &gt; (set Phigh + 5 cmH2O) or monitored Plow &gt; (set Plow +5 cmH2O) for two consecutive breaths.</li> <li>If Tlow is set to &lt; 3 s, the High PEEP alarm is disabled for Plow settings. This reduces the incidence of false positive alarms.</li> </ul>	Check the patient. Check and adjust the ventilator settings, including alarms. Check expiration tubing.

Table 8-2. Alarms and other messages

Alarm	Definition	Action needed
High pressure	<i>High priority.</i> The measured inspiratory pressure > the set High Pressure alarm limit. The ventilator immediately closes the inspiratory valve to stop gas flow to the patient and opens the expiratory valve to reduce pressure to the PEEP/CPAP level. If the pressure reaches 15 cmH <sub>2</sub> O above the High pressure alarm limit for longer than 5 s the ventilator opens the release valve. If the pressure reaches 15 cmH <sub>2</sub> O above the High pressure alarm limit for longer than 5 s ventilator enters the ambient state.	Check the patient. Adjust the Pressure alarm limit. Check the breathing circuit and flow sensor tubes for kinks and occlusions. Provide alternative ventilation once the ventilator enters the ambient state.
Inspiratory volume limitation	<i>Medium priority.</i> The delivered V <sub>t</sub> is > 1.5 times the set V <sub>t</sub> high alarm limit. Pressure is reduced to PEEP level. The APV controls reduce the pressure for the next breath by 3 cmH <sub>2</sub> O. Disabled in noninvasive modes.	Reduce the P <sub>support</sub> setting. Adjust the V <sub>t</sub> high alarm limit.
Invalid option board	<i>Low priority.</i> Installed option board is invalid.	Have the ventilator serviced.
Loss of external power	<i>Low priority, reset when silenced.</i> The HAMILTON-C2 is running on battery power due to loss of its primary power source.	Silence the alarm. Check integrity of connection to primary power source. Check battery status. If you have spare batteries, prepare to swap if necessary. Prepare for possible power loss. Obtain alternative ventilation.
Loss of PEEP	<i>Medium priority</i> <ul style="list-style-type: none"> <li>Pressure during exhalation &lt; (set PEEP/CPAP - 3 cmH<sub>2</sub>O) for more than 10 s</li> <li>Measured end expiratory pressure &lt; (set PEEP/CPAP - 3 cmH<sub>2</sub>O) for two consecutive breaths.</li> </ul>	Check the patient. Check the breathing circuit for leaks. Replace the breathing circuit, if necessary.
Loudspeaker defective	<i>High priority.</i> A loudspeaker malfunction was detected. A technical alarm cannot typically be corrected by the operator. Ventilation continues.	Have the ventilator serviced.

**Table 8-2. Alarms and other messages**

<b>Alarm</b>	<b>Definition</b>	<b>Action needed</b>
Options not found	<i>High priority.</i> Options were not found during startup.	Restart device. If the problem persists, have the ventilator serviced.
Performance limited by high altitude	<i>Medium priority, Low after silence.</i> The airway pressure cannot be reached at the current altitude.  As long as the device remains above the altitude limit, the pressure cannot be reached, and the alarm is active.	Check the patient. Provide alternative ventilation if needed.
Release valve defective	<i>Low priority.</i> During the routine check of the ambient valve during a tightness test, the valve was found to be defective.  The alarm is reset when a tightness test is successfully passed.	If the problem still persists, have the ventilator serviced.
Replace O2 cell	<i>High priority.</i> Communication error, O2 cell defective	Replace O2 cell
Self test failed	<i>High priority.</i> The self test failed during startup. The Start ventilation button is ghosted.  Note that if this error occurs when the device is restarting from a complete power loss, the device enters the ambient state.	Restart device. If the problem persists, have the ventilator serviced.  If the device enters the ambient state, provide alternative ventilation and have the ventilator serviced.
Suctioning maneuver	<i>Low priority.</i> Ventilation suppression is active, and ventilator settings are being maintained, although the ventilator is not delivering breaths.	Resume ventilation when desired by first reconnecting the patient.
Technical state failed	There is a problem with the hardware configuration. Ventilation is impossible.	Have the ventilator serviced.
Unknown part number	There is a problem with the hardware configuration. Ventilation is impossible.	Have the ventilator serviced.

## Appendix A. Specifications

In Appendix A, "Specifications", two sets of tables have been updated:

- Tables A-5 and A-6, "Control setting ranges and resolutions". See "Control settings" on page 9.

- Tables A-12 and A-13, “Adjustable alarm priority, range, defaults, and resolution”. See “Adjustable alarms” on page 10.

### Control settings

Tables A-5 (adult/pediatric) and A-6 (neonatal), “Control setting ranges and resolutions”, have been updated. The following table lists only those settings that have been updated.

**Table A-5. Control settings, ranges and accuracy**

Parameter or Setting (units)	Range		Default settings		Accuracy <sup>1</sup>
	Adult/Ped	Neonatal	Adult/Ped	Neonatal	
ETS <sup>2,3</sup> (%)	5 to 80	5 to 80	25 <i>Noninvasive modes:</i> 35 <sup>3</sup>	25 <i>Noninvasive modes:</i> 35 <sup>3</sup>	
I:E <sup>4</sup>	1:9 to 4:1	1:9 to 4:1	1:4	1:3	--
P-ramp <sup>5</sup> (ms)	0 to 2000  <i>ASV, NIV, NIV-ST, SPONT:</i> max = 200	0 to 600  <i>NIV, NIV-ST, SPONT, nCPAP-PS:</i> max = 200	100	50	±10 ms
Rate <sup>6</sup> (b/min)	<i>(S)CMV+, PCV+:</i> 4 to 80  <i>PSIMV+, NIV-ST:</i> 5 to 80  <i>Other modes:</i> 1 to 80	<i>(S)CMV+, PCV+, PSIMV+, NIV-ST, nCPAP-PS:</i> 15 to 80  <i>PSIMV (non-Intel-lisync):</i> 5 to 80  <i>Other modes:</i> 1 to 150	3.0 to 5.8 <i>IBW:</i> 38  5.9 to 8.0 <i>IBW:</i> 32  8.1 to 20.0 <i>IBW:</i> 25  20.1 to 29.9 <i>IBW:</i> 19  30 to 39 <i>IBW:</i> 17  40 to 59 <i>IBW:</i> 15  60 to 200 <i>IBW:</i> 12	0.2 to 1.25 kg: 60  1.26 to 3.0 kg: 45  3.1 to 5.9 kg: 35 6.0 to 8.9 kg: 30  9.0 to 20.5 kg: 25  21 to 30 kg: 20	Adult: ±1 b/min  Neo: ±1 < 100 ±3 ≥ 100
Tl <sup>4,6,7</sup> (s)	0.1 to 12	0.1 to 12	Based on rate (IBW) and I:E = 1:4	Based on rate (Weight) and I:E = 1:3	±0.01
Tl max <sup>8</sup> (s)	1 to 3	0.25 to 3.0	1.5	1.0 s ≤ 10 kg 1.5 s > 10 kg	± 0.1
VT/kg <sup>9</sup> (ml/kg)	5 to 12	5 to 12	8	5	

1. The stated accuracy includes the tolerance interval for each measurement. See "Accuracy testing" on page 12 for details.
2. Expiratory trigger sensitivity, in % of inspiratory peak flow.
3. When selecting a noninvasive mode, the device uses the ETS value used in the previous mode, if available. If the previous mode did not use ETS, the device sets ETS to 35.
4. In PCV+ and (S)CMV+ modes, mandatory breath timing can be controlled by using a combination of inspiratory time (TI) and rate, or by the I:E ratio; set the method in Configuration. All other modes are controlled by using a combination of inspiratory time (TI) and rate.
5. P-ramp is limited by one-third (1/3) of TI time. Adjustment of TI time can override P-ramp setting.
6. Startup setting derived from body weight setting (neonates), IBW (adults/pediatrics).
7. Inspiratory time; used with Rate to set the breath cycle time.
8. Maximum inspiratory time for spontaneous breaths during noninvasive ventilation.
9. Set in configuration.

### Adjustable alarms

Tables A-12 (adult/pediatric) and A-13 (neonatal), "Adjustable alarm priority, range, defaults, and resolution", have been updated. The following table lists all of the adjustable alarms.

**Table A-12. Adjustable alarm priority, range, defaults, and resolution**

Alarm (units)	Priority	Range		Default setting		Resolution
		Adult/Ped	Neo	Adult/Ped	Neo	
Apnea time <sup>1</sup> (s)	Adult: High Neonatal: Medium	15 to 60	in nCPAP-PS: 5 to 60/ OFF all other modes: 5 to 60	20	15	Adult: 5 s Neonatal: 1 < 15 s 5 ≥ 15
ExpMinVol, low <sup>2</sup> (l/min)	High	in NIV, NIV-ST: OFF/0.1 to 50 other modes: 0.1 to 50	OFF/0.01 to 10	Based on Rate and Vt $0.6 * Rate * Vt$	Based on Rate and Vt $0.6 * Rate * Vt$	Adult: 0.1 < 1 0.5 ≥ 1 1 ≥ 10 Neo: 0.01 < 1 0.1 ≥ 1
ExpMinVol, high <sup>2</sup> (l/min)	High	in NIV, NIV-ST: 0.1 to 50/OFF other modes: 0.1 to 50	0.03 to 10/OFF	Based on Rate and Vt $1.5 * Rate * Vt$	Based on Rate and Vt $1.5 * Rate * Vt$	Adult: 0.1 < 1 0.5 ≥ 1 1 ≥ 10 Neo: 0.01 < 1 0.1 ≥ 1
fTotal, low (b/min)	Medium	0 to 99	0 to 200	0	0	1
fTotal, high (b/min)	Medium	0 to 99	2 to 210	40	70	1

Table A-12. Adjustable alarm priority, range, defaults, and resolution

Alarm (units)	Priority	Range		Default setting		Resolution
		Adult/Ped	Neo	Adult/Ped	Neo	
Oxygen, high <sup>3,4</sup> (%)	High	18 to 105	18 to 105	55	55	1
Oxygen, low <sup>3,4</sup> (%)	High	18 to 97	18 to 97	45	45	1
PetCO <sub>2</sub> , high <sup>5</sup> (mmHg)	Medium	1 to 100	1 to 100	60	60	1
(kPa)		1 to 13.2	1 to 13.2	8	8	0.1
PetCO <sub>2</sub> , low <sup>5</sup> (mmHg)	Medium	OFF/0 to 100	OFF/0 to 100	30	30	1
(kPa)		OFF/0 to 13.2	OFF/0 to 13.2	4	4	0.1
Pressure, high (Pmax) (cmH <sub>2</sub> O)	High	15 to 70	15 to 70	40	40	1
Pressure, low (cmH <sub>2</sub> O)	High	4 to 60	<i>nCPAP-PS:</i> 2 to 60 <i>other modes:</i> 4 to 60	PEEP	<i>nCPAP-PS:</i> 2 <i>other modes:</i> PEEP	1
Pressure limitation (cmH <sub>2</sub> O)	Medium; Low after silence	5 to 60	5 to 60	Pmax - 10	Pmax - 10	1

**Table A-12. Adjustable alarm priority, range, defaults, and resolution**

Alarm (units)	Priority	Range		Default setting		Resolution
		Adult/Ped	Neo	Adult/Ped	Neo	
Vt, low <sup>6</sup> (ml)	Medium	OFF <sup>7</sup> /10 to 3000	OFF <sup>7</sup> /0.1 to 300	Based on IBW $0.5 * Vt$	Based on Weight $0.5 * Vt$	<i>Adult:</i> OFF 5 < 100 ml, 10 ≥ 100 and < 500, 50 ≥ 500 <i>Neo:</i> OFF 0.1 < 10, 1 ≥ 10 and < 100, 5 ≥ 100
Vt, high <sup>6</sup> (ml)	Medium	10 to 3000/ OFF	0.1 to 300/OFF	Based on IBW $1.5 * Vt$	Based on Weight $1.5 * Vt$	<i>Adult:</i> OFF 5 < 100 ml, 10 ≥ 100 and < 500, 50 ≥ 500 <i>Neo:</i> OFF 0.1 < 10, 1 ≥ 10 and < 100, 5 ≥ 100

1. The default setting is configurable.
2. Startup setting derived from body weight setting (neonates), IBW (adults/pediatrics).
3. Active only when O2 monitoring (O2 sensor) is enabled.
4. The high and low oxygen alarm limits are automatically set in relation to the current oxygen setting as follows: O2 setting + 5 (Oxygen high limit) and O2 setting - 5 (Oxygen low setting). For example, if the Oxygen setting is 70%, the Oxygen high limit is set to 75 and the low limit is set to 65.
5. Requires CO2 option.
6. In ASV mode, this alarm only applies for spontaneous breaths.
7. OFF available in noninvasive and neonatal modes.

## Accuracy testing

The ventilator's parameter and measurement accuracy is tested using an IMT FlowAnalyser™. The tolerance intervals for the data generated by the FlowAnalyser are as specified below, and are included in the accuracy information provided in this manual.

**Table A-13. Tolerance intervals for accuracy testing**

Parameter type	Tolerance interval of measurement
Volume	≤ 50 ml: ±1% > 50 ml: ±1.75%
Pressure	±0.75% or ±0.1 cmH2O, whichever is greater
Flow	±1.75% or ±0.5 l/min, whichever is greater

**Table A-13. Tolerance intervals for accuracy testing**

Parameter type	Tolerance interval of measurement
O <sub>2</sub>	±1 %

## Appendix J. Configuration

Two areas in Configuration have been changed:

- The gas source (HPO or LPO) is no longer set in Configuration. See “Selecting the gas source (HPO, LPO)” on page 13.
- You can now set a minimum alarm loudness for the device. See “Setting the minimum alarm loudness” on page 13.

### Selecting the gas source (HPO, LPO)

The method for selecting the gas source to use, HPO or LPO, has changed. The controls have been moved from Configuration to the Utilities window, and can be changed in Standby.

#### To select the oxygen source

- 1 In Standby mode, touch the **Utilities** button. By default, the Gas source window is displayed.
- 2 Touch the appropriate button for the desired oxygen source.
  - Select **HPO** for high-pressure oxygen (the default)
  - Select **LPO** for low-pressure oxygen

The ventilator always resets to HPO mode when restarted.
- 3 Close the Utilities window.

### Setting the minimum alarm loudness

You can set a minimum alarm loudness (volume) setting for the device. Once set, the device operator cannot set the alarm volume below the value set here in Configuration.

#### To set the minimum alarm loudness (volume)

- 1 In Configuration, open the General -> More window.
- 2 Touch the **Min. Loudness** button and choose the minimum alarm volume to allow on the device. By default, it is set to 1.
- 3 Continue setting configuration options or exit Configuration mode.

The setting is applied to the device. Note that if the new minimum is greater than the currently set alarm volume, the alarm volume is reset to the new minimum level. To verify the setting, check the **Loudness** value in the System -> Settings window.

Hamilton Medical AG  
Via Crusch 8, 7402 Bonaduz,  
Switzerland  
☎ +41 58 610 10 20  
info@hamilton-medical.com  
www.hamilton-medical.com

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