
<b>Instructions for use</b>	Breathing circuit set, coaxial .....	4
<b>Gebrauchsanweisung</b>	Beatmungsschlauchset, koaxial .....	8
<b>Instructions d'utilisation</b>	Kit de circuit respiratoire, coaxial .....	13
<b>Instrucciones de uso</b>	Equipo respiratorio coaxial .....	18
<b>Istruzioni per l'uso</b>	Set circuito paziente, coassiale .....	23
<b>Instruções de uso</b>	Kit de respiração, coaxial .....	28
<b>使用说明</b>	呼吸装置, 同轴 .....	33
<b>使用説明書</b>	同軸呼吸回路セット .....	37
<b>Инструкции по эксплуатации</b>	Комплект дыхательного контура, коаксиальный .....	41
<b>Kullanım talimatları</b>	Solunum devresi seti, koaksiyel .....	46

**REF** 260087, 260094, 260145, 260144, 260128, 260127, 260167, 260168

CE 0197


ERC

**HAMILTON**  
**MEDICAL**

- |   |   |
|---|---|
| <p><b>1</b> To patient<br/>Zum Patienten<br/>Port inspiratoire<br/>Hacia el paciente<br/>Al paziente<br/>Para o paciente<br/><b>至病人</b><br/>吸気側接続口<br/><b>Порт «к пациенту»</b><br/>Hastaya</p> <p><b>2</b> From patient<br/>Vom Patienten<br/>Port expiratoire<br/>Desde el paciente<br/>Dal paziente<br/>Do paciente<br/><b>自病人</b><br/>呼気側接続口<br/><b>Порт «от пациента»</b><br/>Hastadan</p> <p><b>3</b> Coaxial inspiratory/expiratory limb<br/>Koaxialer Inspirations-/Expirationsschenkel<br/>Branche coaxiale inspiratoire / expiratoire<br/>Rama inspiratoria/espíratória coaxial<br/>Branca inspiratoria/espíratória coassiale<br/>Alça inspiratória/espíratória coaxial<br/><b>同軸吸気肢/呼気肢</b><br/><b>同軸吸気/呼気回路</b><br/><b>Коаксиальный патрубок вдоха/выдоха</b><br/>Koaksiyel inspiratuvar/ekspiratuvar parça</p> | <p><b>4</b> Flow sensor (FS)<br/>Flow-Sensor (FS)<br/>Capteur de débit (FS)<br/>Sensor de flujo (FS)<br/>Sensore di flusso (FS)<br/>Sensor fluxo (FS)<br/><b>流量传感器 (FS)</b><br/>フローセンサ (FS)<br/><b>Датчик потока (FS)</b><br/>Akış sensörü (FS)</p> <p><b>5</b> Expiratory valve set (EV)*<br/>Expirationsventil (EV)*<br/>Système de valves expiratoires (EV)*<br/>Conjunto de válvula espiratoria (EV)*<br/>Set valvola espiratoria (EV)*<br/>Conjunto de válvula expiratória (EV)*<br/><b>呼气阀套件 (EV)*</b><br/><b>呼気弁セット (EV) *</b><br/><b>Комплект клапана выдоха (EV)*</b><br/>Ekspiratuvar valf seti (EV)*</p> |
|---|---|

\* Not included in all breathing circuit sets | Nicht bei allen Beatmungsschlauchsets im Lieferumfang enthalten | Non inclus dans tous les kits de circuit respiratoire | No incluido en todos los equipos respiratorios | Non incluso in tutti i set circuito paziente | Não incluído em todos os kits de respiração | 并非所有呼吸装置中都包括 | すべての呼吸回路セットには含まれていません | **Включен не во все комплекты дыхательных контуров.** | Tüm solunum devresi setlerinde yer almaz

REF

260087 (1.8 m + FS), 260094 (2.4 m + FS), 260145 (3.0 m + FS), 260144 (4.8 m + FS)

REF

260128 (1.8 m + FS + EV), 260127 (2.4 m + FS + EV), 260167 (3.0 m + FS + EV), 260168 (4.8 m + FS + EV)

## Conventions

### ⚠ WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

### ⚠ CAUTION

*Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.*

### NOTICE

Emphasizes information of particular importance.

## 1 Overview

The coaxial breathing circuit set with flow sensor (PN 260087, 260094, 260145, 260144) as well as the coaxial breathing circuit set with flow sensor and expiratory valve (PN 260128, 260127, 260167, 260168) can be used with HAMILTON-C1/T1/MR1 ventilators.

### 1.1 Intended use

#### PN 260087, 260094

The ventilator breathing circuits adult/pediatric with PN 260087 and 260094 are intended to be used to connect the HAMILTON-C1/T1/MR1 ventilators to the patient tube or respiratory mask during ventilation.

The Hamilton Medical flow sensor is a device to measure patient air flow with Hamilton Medical ventilators.

All components of the breathing set are intended for use by qualified, trained healthcare professionals.

#### PN 260145, 260144

The ventilator breathing circuits adult/pediatric with PN 260145 and 260144 are intended to be used to connect the HAMILTON-C1/T1/MR1 ventilators to the patient tube or respiratory mask during ventilation only in environments where shorter breathing circuits are not applicable (for example, with MR/CT scanners).

The Hamilton Medical flow sensor is a device to measure patient air flow with Hamilton Medical ventilators.

All components of the breathing set are intended for use by qualified, trained healthcare professionals.

#### PN 260128, 260127

The ventilator breathing circuits adult/pediatric with PN 260128 and 260127 are intended to be used to connect the HAMILTON-C1/T1/MR1 ventilators to the patient tube or respiratory mask during ventilation.

The Hamilton Medical flow sensor is a device to measure patient air flow with Hamilton Medical ventilators.

The Hamilton Medical single-use adult/pediatric expiratory valve set for HAMILTON-C1/T1/MR1 is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended for use with equipment other than Hamilton Medical ventilators.

All components of the breathing set are intended for use by qualified, trained healthcare professionals.

#### **PN 260167, 260168**

The ventilator breathing circuits adult/pediatric with PN 260167 and 260168 are intended to be used to connect the HAMILTON-C1/T1/MR1 ventilators to the patient tube or respiratory mask during ventilation only in environments where shorter breathing circuits are not applicable (for example, with MR/CT scanners).

The Hamilton Medical flow sensor is a device to measure patient air flow with Hamilton Medical ventilators.

The Hamilton Medical single-use adult/pediatric expiratory valve set for HAMILTON-C1/T1/MR1 is an accessory for the exhalation gas path. Its main function is to control expiratory gas flow and maintain user-set pressure levels. The Hamilton Medical single-use expiratory valve set is not intended

for use with equipment other than Hamilton Medical ventilators.

All components of the breathing set are intended for use by qualified, trained healthcare professionals.

## **2 Safety information**

### **⚠ WARNING**

- Reprocessing of Hamilton Medical single-use products can affect the product properties and may cause injury to the patient. For example, a change to the surface structure during reprocessing may lead to a change in the tear strength or cause actual cracking. Furthermore, an altered surface structure may result in, for example, a microbial aggregation of spores, allergens, and pyrogens, or cause an increase in the number of particles released as a result of chemical changes in the material properties.
- Hamilton Medical does not assume any warranty for the proper function of single-use items if they are reprocessed and reused by the user. Any liability on the part of Hamilton Medical for any damages resulting from the use of reprocessed or reused single-use parts, or use of any non-Hamilton Medical accessories and consumables, is excluded.
- Follow hospital infection control procedures for all tests and use.

- Remove the safety cover(s) before use.
- Discard the product if it fails two subsequent preoperational checks or if there is any sign of damage.
- To prevent patient or ventilator contamination, always use a bacteria filter or HMEF between the patient and the inspiratory port. If no bacteria filter is used, the exhaled gas can contaminate the ventilator.
- Adding attachments or other parts/assemblies to a breathing system can change the pressure gradient across the ventilator, which can adversely affect ventilator performance.
- Connect only a single breathing circuit at a time. Do not connect any circuit as an extension to another one. Doing so can seriously harm or kill the patient.
- Route the breathing circuit without tension and without any kinks from the ventilator or humidifier to the patient, and protect against unintended influences.
- Mount the expiratory valve set outside of the MR environment only! Otherwise the MRI scanner and the valve membrane may be severely damaged.
- When used with pediatric patients, the compliance of the 4.8 meter breathing circuits (PN 260144, 260168) may lead to increased pressures or tidal volumes, depending on the ventilation mode.

## ⚠ CAUTION

- (USA only) Federal law restricts this device to sale by or on the order of a physician.
- Handle used products as contaminated goods according to local laws and regulations or hospital internal procedures.
- Always carry out all preoperational checks before connecting the patient to prevent possible patient injury. Make sure another source of ventilator support is available.

## NOTICE

- Before use, read the appropriate ventilator *Operator's Manual*.
- Replace the product in accordance with hospital infection control procedures, or depending on the patient's secretions and nebulization of medication. Consider using a filter at the expiratory port during nebulization.
- Any bacteria filter, HMEF, or additional accessories in the expiratory limb may substantially increase flow resistance and impair ventilation.
- When adding parts to the Hamilton Medical breathing circuit configurations, do not exceed the inspiratory and expiratory resistance values of the ventilator breathing system as specified in the ventilator *Operator's Manual*.

- Certain parts of the breathing circuit may be visible on the MRI image. Keep the breathing circuit out of the scanned area as much as possible.
- To prevent inaccurate flow sensor readings, make sure the flow sensor is correctly connected: The flow sensor tubes must not be kinked, and must be secured with the included clamp.
- Ensure the expiratory valve set is in its intended position in the ventilator.
- Ensure you select the correct expiratory valve (adult/pediatric) for your patient. If the expiratory valve type does not match the selected patient group on the ventilator, the **Wrong expiratory valve** alarm is generated.
- Ensure the membrane is aligned with the expiratory valve body and the metal faces up.

## 3 Connecting and calibrating components

### To connect the breathing circuit


- ▶ Connect the breathing circuit as shown on page 3, and ensure that the limbs are connected to the appropriate ports on the ventilator.

### To connect the flow sensor

- 1 Remove the safety cover.
- 2 Connect the flow sensor to the breathing circuit. Ensure the end marked **PATIENT** faces the patient.
- 3 Connect the blue and clear tubes to the same-colored flow sensor connectors on the ventilator. Blue to blue, clear to white/silver.
- 4 Calibrate the flow sensor.

### To install the expiratory valve

- 1 Remove the safety cover.
- 2 Ensure the membrane is properly aligned with the expiratory valve housing and the metal plate faces up (**A**).
- 3 Position the expiratory valve set in the expiratory port (**B**) and twist the locking ring clockwise until it locks into place (**C**).

## 4 Preoperational check: Leak test and calibration

Perform all of the required preoperational checks before patient use, including a Leak test and necessary calibrations. For detailed instructions, see your ventilator *Operator's Manual*.

## 5 Delivery

The components are delivered clean and ready for clinical use.

## 6 Disposal

A used component must be handled as contaminated. Follow all local, state, and federal regulations with respect to waste management and environmental protection when disposing of used parts.

## 7 Environmental requirements and storage conditions

### Operating temperatures and storage conditions

Operating temperature	-15°C to 50°C
Storage temperature	-20°C to 60°C
Transport	-20°C to 60°C (max 4 weeks)
Storage and operating humidity	5% to 95% relative humidity, noncondensing

## 8 Technical specifications

PN 260087, 260094, 260145, 260144, 260128, 260127, 260167, 260168	
Compliance for 1.80 m for 2.40 m for 3.0 m for 4.80 m	1.7 ml/cmH <sub>2</sub> O at 60 cmH <sub>2</sub> O 2.4 ml/cmH <sub>2</sub> O at 60 cmH <sub>2</sub> O 2.7 ml/cmH <sub>2</sub> O at 60 cmH <sub>2</sub> O 4.7 ml/cmH <sub>2</sub> O at 60 cmH <sub>2</sub> O
Leakage	< 50 ml/min
Conforms to	EN 12342 / ISO 5367:2014
Tube length	1.8 m / 2.4 m / 3.0 m / 4.8 m
Conical tube connectors	Conforms to EN ISO 5356-1
Smallest internal diameter	13 mm
Max. operating duration	28 days
Conforms to	ISO 5367:2014
Maximum working pressure	200 cmH <sub>2</sub> O

Resistance to flow at	15 l/min	30 l/min
<b>260087, 260128 (1.80 m)</b>		
Inspiratory (cmH <sub>2</sub> O)	0.01	0.02
Expiratory (cmH <sub>2</sub> O)	0.02	0.04
<b>260094, 260127 (2.4 m)</b>		
Inspiratory (cmH <sub>2</sub> O)	0.02	0.03
Expiratory (cmH <sub>2</sub> O)	0.03	0.04
<b>260145, 260167 (3.0 m)</b>		
Inspiratory (cmH <sub>2</sub> O)	0.03	0.04
Expiratory (cmH <sub>2</sub> O)	0.03	0.05
<b>260144, 260168 (4.8 m)</b>		
Inspiratory (cmH <sub>2</sub> O)	0.03	0.06
Expiratory (cmH <sub>2</sub> O)	0.03	0.06

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Specifications are subject to change without notice. For all proprietary and third-party trademarks used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). © 2021 Hamilton Medical AG. All rights reserved.
Hamilton Medical AG
medin Medical Innovations GmbH