Reprocessing Guide

Flow sensor

REF 950185

624421/05 | 2017-09-30

CE 0197

HAMILTON
MEDICAL

Intelligent Ventilation since 1983

Table of Contents

English	3
Deutsch	13
Français.....	23
Español.....	33
Português	43
Italiano.....	53
中文.....	63
日本語	71
Русский	81
Türkçe	93

Reprocessing guide per EN ISO 17664 (ENGLISH)

© 2017 Hamilton Medical AG. All rights reserved. Printed in Switzerland. No part of this publication may be reproduced or stored in a database or retrieval system, nor transmitted, in any form or by any means, electronic, mechanical, by photocopying, recording, or otherwise, without the prior written permission of Hamilton Medical.

This manual may be revised or replaced by Hamilton Medical at any time and without notice. Ensure that you have the most current applicable version of this manual; if in doubt, contact Hamilton Medical AG Marketing Department. While the information set forth is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this manual shall limit or restrict in any way Hamilton Medical's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated and serviced by trained professionals only.

Definitions

WARNING

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

CAUTION

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

NOTE

Emphasizes information of particular importance.

Table of Contents

Definitions	3
Reprocessing after use and preparation	6
Manual cleaning and disinfection	6
Mechanical cleaning and disinfection	7
Packaging	9
Visual test	9
Sterilization	10
Testing before use	10
Flow sensor life span	10
Autoclaved and packaged flow sensor: sterility and storage conditions	10
Disposal	11

NOTE

- The autoclavable flow sensor consists of the following materials: Polysulfone (housing), stainless steel (flap), and silicone (measuring tubes). All materials used are heat resistant up to 140°C (284°F).
- The tube attachment strap consists of the following materials: Silicone 50 shore A
- Reprocess the tube attachment strap and the calibration adapter the same way the flow sensor is reprocessed.

WARNING

- **Before reusing the flow sensor, it must be cleaned, disinfected, and sterilized.**
- **Hamilton Medical cannot be held responsible for the correct functioning of flow sensors that are not reprocessed according to these instructions.**
- **Ensure that only processes that have been specifically validated for the product or device are used, and that the validated parameters are used with every cycle.**
- **The flow sensor must be handled as a contaminated item. Follow all local, state, and federal regulations with respect to environmental protection when disposing of used flow sensors.**
- **Follow hospital infection control procedures, as well as local laws and regulations. This applies in particular to the various regulations regarding an effective deactivation of prions.**
- **Visually inspect the flow sensor body, tubing, and internal flap.**
- **Discard the flow sensor if there is any sign of damage.**

CAUTION

- The autoclavable flow sensor has a limited life span.
 - If reprocessing flow sensors manually, take care that no damage is caused by the use of hard brushes, scouring agents, or by the exertion of too much force. Do not insert brushes into the flow sensor, as this will damage the flap.
 - Note that the use of rinse aids will reduce the life span of the flow sensor. Rinse aids, such as Neodisher Mediklar, lead to early failure and cracks in the plastic flow sensor body.
 - Calibrate the flow sensor prior to patient use. If calibration fails, it may be repeated once. The flow sensor must be discarded if the calibration fails the second time.
 - Reprocess the flow sensor's calibration connector the same way the flow sensor is reprocessed.
-

Reprocessing after use and preparation

Used flow sensors may be contaminated. The flow sensor must be reprocessed immediately after use. The reprocessing process cleans, disinfects, and sterilizes the flow sensor.

The flow sensor must be cleaned and disinfected prior to sterilization. Due to considerably higher efficiency, Hamilton Medical recommends that you always use a mechanical process (a disinfectant) for cleaning and disinfection.

Ultrasound cleaning

When large impurities or matter are present in the flow sensor body or tubes, you must clean the components in an ultrasound bath prior to normal cleaning.

- 1 Soak the flow sensor in the ultrasound bath for five minutes.
The water in the ultrasound bath must not be hotter than 40°C (104°F).
You can optionally add a cleaning agent (Neodisher Mediclean forte 0.5 - 2%).
- 2 Remove macroscopic impurities of the flow sensor body and tubes by rinsing or wiping.
You can add a cleaning agent (Neodisher Mediclean forte 0.5 - 2%) to the rinse water.
You must not use any hard tools or hard brushes to remove resilient impurities.
- 3 Rinse the flow sensor, including both tubes, at least five times intensively, or according to the validated cleaning plan, in 30 ml freshly distilled or deionized water.

Manual cleaning and disinfection

NOTE

As mechanical disinfection is more effective and consistent, manual cleaning and disinfection is only permitted when no mechanical process is available.

When selecting the cleaning and disinfection agent, consider whether the agents in question are suitable for the flow sensor, and are compatible with each other. Make sure the disinfection agents' effects are proven and the chemicals are compatible with the materials of the flow sensor. In addition, instructions for cleaning with the selected agents are required.

Follow the chemical concentrations and soak times as stated in the manufacturer's Instructions for use. Only use freshly made solutions. The disinfection solution must not foam. Use only sterile water or water with a low microorganism count for all cleaning steps. Make sure that the water contains no endotoxins and contains minimal particle pollutants.

Manual cleaning

- 1 Clean the components in an ultrasound bath prior to normal cleaning, as described in "Ultrasound cleaning" on page 6.
- 2 Submerge the flow sensor in the cleaning solution and let it soak for the time defined by the manufacturer of the disinfection or cleaning agent (Cidezyme). Make sure that all parts of the flow sensor are fully submerged in the solution.
- 3 Rinse both tubes at the beginning and the end of the soak time at least five times with the cleaning agent.
- 4 Remove matter and larger exterior impurities by carefully scrubbing the flow sensor with a soft brush or soft towel.
- 5 Rinse the flow sensor, including both tubes, at least five times intensively, or according to the validated cleaning plan, in freshly distilled or deionized water. Repeat the cleaning process if the last cleaning solution was not clear or there are still visible impurities on the flow sensor.

Manual disinfection

- 1 Submerge the flow sensor in the disinfection solution and let it soak for the time defined by the manufacturer of the disinfection agent (CIDEX OPA). Make sure that all parts of the flow sensor are fully submerged in the solution.
- 2 Rinse both tubes at the beginning and the end of the soak time at least five times, or in accordance with the validated disinfection plan, with the disinfection solution.
- 3 Rinse the flow sensor, including both tubes, at least five times intensively, or according to the validated cleaning plan, in freshly distilled or deionized water. Repeat the cleaning process if the last cleaning solution was not clear or there are still visible impurities on the flow sensor.
- 4 Dry the flow sensor with filtered, oil-free compressed air. Do not apply the compressed air directly to the flow sensor's flap as it might cause damage.
- 5 Immediately package the flow sensor using appropriate packaging material.


Mechanical cleaning and disinfection

The products must be reprocessed in such a manner that hygienic and safe reuse can be assured. Cleaning should only be carried out in a cleaning and disinfection device that complies with ISO 15883 and has been proven to be effective.

The disinfectant must have cleaning connectors (injectors) to connect to lumina (e.g., cannulae, handles) that have sufficient and reproducible cleaning pressure (confirmed in the framework of a specific validation).

- 1 Clean the components in an ultrasound bath prior to normal cleaning, as described in "Ultrasound cleaning" on page 6.

- 2 Connect each of the flow sensor tubes to an injector. Hamilton Medical recommends the E453 (4 X 110 mm) injector with enclosed retaining clamp.

- 3 Insert the flow sensor tubes into the device without bending or kinking, while leaving the flow sensor exposed to the cleaning and disinfection agents. Make sure the tubes do not sag. In addition, when affixing the flow sensors, make sure that they do not impede the effectiveness of the cleaning and disinfection. Flow sensors must not touch each other.
- 4 Place the tube attachment strap in an autoclave wire basket as shown above. Make sure you lock the basket.
- 5 When complete, package the flow sensor and sensors using appropriate packaging material.

The following program parameters must be met for successful mechanical cleaning/disinfection:

Pre-rinse:	one cycle using cold water for 1 min
Cleaning:	one cycle at 55°C (131°F) for 5 min
Optional neutralization:	one cycle using cold water for 1 min
Rinsing:	one cycle using cold water for 1 min

Thermic disinfection:	one cycle at 83°C (181.4°F) for 10 min
Drying:	100°C (212°F) for 10 min and 95°C (203°F) for 30 min

Follow the instructions provided in the operator's manual of the cleaning and disinfection device. To ensure safe internal cleaning, the flow sensor must be connected to the corresponding injectors. The flow sensor must not disconnect from the injector during reprocessing.

Flow sensors that disconnect during reprocessing must be processed again. After the cleaning process is complete, check that the flow sensor is completely dry and undamaged. Damaged flow sensors must be replaced.

Recommended equipment

Hamilton Medical recommends the DES-VAR-TD-Anaesthesia program, among others in the Miele PG8536 disinfectant, together with the E436/6 injector tray, and the E453 Miele injectors.

Suitable cleaning agents:

Manufacturer	Product	Concentration
Dr. Weigert	Mediclean forte®	1.00%

Suitable neutralizer

Manufacturer	Product	Concentration
Dr. Weigert	Neodisher Z®	0.1%

CAUTION

Using a rinse agent will cause premature damage and reduce product life span.

Packaging

Make sure that the flow sensors are not moist during packaging.

The packaging must conform to ISO 11607 and be suitable for vapor sterilization (heat resistance up to 141°C (285.8°F)) and be sufficiently permeable to vapor.

Only use packaging suitable for sterilization.

Visual test

After each cleaning and disinfection cycle, the flow sensor must be macroscopically clean, that is, free of visible residual matter and other impurities. If it is not, the entire cleaning and disinfection process must be repeated.

Visually check for external damage, such as cracks, broken or deformed parts, or discoloration.

Sterilization

Sterilize the flow sensor after cleaning and disinfection before use. Use one of the following methods:

- 134°C (273.2°F) with or without prevacuum, with an exposure time of minimum 3 min and maximum 18 min
- 132.2°C (270°F) with or without prevacuum, with an exposure time of minimum 3 min and maximum 18 min
- 121°C (249.8°F) with or without prevacuum, with an exposure time of minimum 30 min

Other sterilization methods, such as hot-air, ethylene oxide, formaldehyde, radiation, and low-temperature plasma sterilization are outside Hamilton Medical's responsibility.

Testing before use

WARNING

Defective flow sensors or flow sensors with faulty calibration must not be used.

Calibrate the flow sensor before use on the patient. To calibrate, follow the instructions provided in the appropriate ventilator operator's manual. If calibration fails, it may be repeated once. The flow sensor must be discarded if the calibration fails the second time.

Before using the flow sensor, carry out a visual test and a tightness test, as described in the ventilator operator's manual and in the *Autoclavable Flow Sensor Instructions for Use* (PN 624610). Discard defective flow sensors.

Flow sensor life span

Discard damaged or deformed flow sensors, as well as flow sensors that fail two subsequent calibrations.

Autoclaved and packaged flow sensor: sterility and storage conditions

The sterility of an autoclaved and packaged flow sensor depends on how long the packaging can keep the flow sensor sterile. Follow the packaging manufacturer's specifications. At a minimum, the flow sensor must be autoclaved every two years. Storage is subject to the Hamilton Medical guidelines as specified in the *Autoclavable Flow Sensor Instructions for Use* (PN 624610).

Disposal

The flow sensor must be handled as a contaminated item. Follow all local, state, and federal regulations with respect to environmental protection when disposing of used flow sensors.

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